



All conversations.  
All mobile phones.  
All the freedom  
of rechargeable.

Moxi All




Driven by Tempus

unitron Hearing matters

A Sonova brand





Moxi All

Product shown at actual size



#### Made for all phones

Moxi All hearing instruments  
connect directly to all phones\*  
– not just iPhones.®



#### Hands-free

There's no need to touch  
the phone. The patient  
clearly hears what the caller  
is saying, while microphones  
transmit their voice to their  
mobile phone.




#### No streamer

Moxi All acts like a Bluetooth®  
headset, connecting to mobile  
phones up to 10m away. Patients  
experience amazing sound  
quality and more enjoyable  
conversations.

\*Smartphones and traditional mobile phones with a compatible Bluetooth Hands-Free Profile  
The Bluetooth word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.

# “You’ll never guess where I’m calling from...”

Because life’s best conversations  
happen on mobile phones too





Moxi All R

# The rechargeable revolution continues



100%

say patients like an overnight charge that lasts all day\*

With the arrival of Moxi All R, patients have the freedom to stay in the conversation all day on a single overnight charge. And patients never have to worry if they forget to charge, because they have the flexibility to swap in traditional batteries at any time.



**Power to your patients!**

Fully charged in less than 7 hours

On average 16 hours of use including 2 hours of phone calls and 5 hours of TV


100% recyclable battery

Available in 6 colors

\*100% of hearing healthcare professionals agree patients like that the hearing instruments are fully charged overnight and last for an entire day.  
Source: Anovum (2017) Unitron Launch Monitor

# Solve their biggest listening challenges. **Everywhere.**

Binaural processing SoundCore features enable patients to experience fantastic speech understanding and amazing spatial awareness thanks to the technology in the Tempus™ platform.



**SoundCore™**  
Best in conversations

Experience automatic, seamless adjustments to **hear their best** across all listening situations, courtesy of SoundNav.

**Understand speech and know where it's coming from**, even in the most challenging conversations in noise and crowds, thanks to binaural SpeechPro.

**Experience clearer speech and more natural listening** everywhere with Sound Conductor.

**Enjoy a more natural, realistic sound experience** with binaural Spatial Awareness.

## **Latest technology keeps them connected**

Our high-end binaural features are fully functioning in Moxi All thanks to the SWORD™ 2.4 GHz chip. It also powers the direct connectivity to mobile devices with Bluetooth and to the TV Connector with AirStream™ Technology.

## **Easy programming**

Programming Moxi All hearing instruments with Noahlink Wireless means patients don't need to wear an extra device.

# 89%

say their patients hear better in conversations with Tempus.\*



\*89% of hearing healthcare professionals said their patients hear better in conversations with Tempus than with any other hearing instrument they have tried. Source: Anovum (2017) Unitron Launch Monitor



# Take personalization to the next level with FLEX

Today's informed healthcare consumers want to feel more connected, confident and in control. We want to make them happy. It's a single-minded focus that's driving the creation of our FLEX™ ecosystem. Although it's named after our flagship FLEX solutions, the FLEX ecosystem is so much more than a portfolio of innovations. It's a world view centered around giving consumers the easy, inclusive experience they crave.

**Because when they feel empowered, we all win.**



87%

of patients consider a trial beneficial<sup>1</sup>

99%

say **FLEX:TRIAL™** improves the total experience<sup>2</sup>

85%

of patients say an upgrade feature would make purchase decisions easier<sup>1</sup>

<sup>1</sup> Source: Consumer Research, Sonova Market Insights (2017)

<sup>2</sup> Source: Anovum (2017) Unitron Launch Monitor

<sup>3</sup> The Remote Plus app for Moxi All works on Apple smartphones with iOS 10.2 or newer and Android smartphones with version 6 or newer. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Google Play and the Google Play logo are trademarks of Google Inc.

## Give them confidence


Patients can take home hearing instruments to try before they buy with **FLEX:TRIAL**. And they have peace of mind that their hearing instruments can meet all their future needs thanks to the flexibility of **FLEX:UPGRADE™**

## Know before it's too late

Get an early warning email if your patient is making negative ratings and quickly see the results online with myInsights.

## Unlock patient insights

Your patients get a much clearer picture of their needs from the combination of Patient Ratings and listening lifestyle evidence gathered by Log It All. You can use this evidence to increase your patients' confidence in your recommendations.



**Remote Plus app**

## Give them a voice

It's easy for patients to share their impressions right in the moment with **Patient Ratings**. They simply download the Remote Plus app<sup>3</sup> to their smartphone and they're ready to make ratings everywhere they go, with no streamer to wear, pair or charge.







# Designed to be easy

Moxi All hearing instruments are designed to put patient needs first, delivering ultimate comfort, ease of use and freedom.



**Easily connect them to TV**  
Our TV Connector accessory automatically switches Moxi All hearing instruments to the TV program, and frees patients from wearing a streamer. Patients instantly enjoy high-quality stereo sound, wirelessly through their hearing instruments, at a comfortable volume for everyone. And it works with laptops, tablets and stereos too.



## Moxi All

Something remarkable for everyone

01 Beige | P2 Amber | S2 Amber suede | Q9 Cinnamon | S6 Sand storm | S3 Espresso boost | S7 Teal blast | P6 Platinum | S5 Pewter shine  
P7 Pewter | P4 Espresso | P8 Charcoal

Products shown at actual size  
\*Colors available for Moxi All R

# Moxi All at a glance

SoundCore		T Pro	T 800	T 700	T 600	T 500
<b>SoundNav</b> Automatic program identifies and classifies sounds across environments, with a focus on conversations.	Music	•	•			
	Noise	•	•	•		
	Conversation in a crowd	•	•	•		
	Conversation in a small group	•	•	•		
	Conversation in quiet	•	•	•	•	
	Conversation in noise	•	•	•	•	
	Quiet	•	•	•	•	
<b>SpeechPro</b> Provides optimal speech understanding and localization, even in the most challenging environments.	Total environments	7	7	6	3	AutoMic
	SpeechPro	•	SpeechZone 2	SpeechZone		
	Speech Locator	•	•	•		
	Speech Focus	•				
<b>Sound Conductor</b> Dynamically balances features to provide speech understanding, comfort and natural sound.	Dynamic Spatial Awareness	•				
	Speech enhancement	•	•	•	•	•
	Noise reduction	•	•	•	•	•
	Adaptive directionality	Multiband	Multiband	Multiband	Multiband	•
<b>Spatial Awareness</b> Makes it easier to localize sound, providing a more natural, realistic listening experience.	Spatial Awareness	Dynamic	Personalized	•		
	Pinna Effect	•	•	•	•	•
<b>Sound Stabilization</b>						
	AntiShock 2	•	•	•	•	•
	Wind control	•	•	•	•	•
	Feedback manager	•	•	•	•	•
	Natural Sound Balance	•	•	•	•	•
<b>Experience innovations</b>						
<b>Patient insights</b>	Log It All	•	•	•	•	•
	Patient Ratings	•	•	•	•	•
	Data logging	•	•	•	•	•
<b>FLEX</b>	FLEX:TRIAL	•	•	•	•	•
	FLEX:UPGRADE		•	•	•	•
<b>Convenience</b>						
	Direct connectivity	•	•	•	•	•
	Rechargeable option	•	•	•	•	•
	DuoLink	•	•	•	•	•
<b>Fitting</b>						
	Automatic Adaptation Manager	•	•	•	•	•
	MyMusic	Binaural	Binaural	•	•	•
	Frequency compression	•	•	•	•	•
	Tinnitus masker	•	•	•	•	•
	IntelliVent	•	•	•	•	•
	Streaming programs	•	•	•	•	•
	Manual programs	•	•	•	•	•
	NAL-NL2/NL1 and DSLv5	•	•	•	•	•
	Fitting channels	20	20	16	10	6

Moxi All is rated IP 68 for protection against moisture and debris.

Technology level availability varies by market.

At Unitron, we care deeply about people with hearing loss. We work closely with hearing healthcare professionals to provide hearing solutions that improve lives in meaningful ways. Because hearing matters.

© 2018 Unitron. All rights reserved.

2018-02 028-6283-02

**unitron.com**

**sonova**  
HEAR THE WORLD